



Release notes LinnéFiler 4.2

Version: 2015-08-27

Author: Camilla Lyttkens, Fertsoft AB

Purpose

These release notes concerns the changes from the previous version that are implemented in version 4.2 of LinnéFiler and all subsystems that are included in the product, where Fertsoft AB is responsible for development and configuration.

System requirements version 4.2 and later

This release of LinnéFiler requires FileMaker Pro version 13 or later as client and FileMaker Server 13 or later as server application. For more information about hardware and operation system requirements, see www.filemaker.com

SYSTEM

Multiple windows

Extended support for multiple windows.

Speed

The overall speed of the system is now increased due to optimization of time-consuming processes.

ivfServices

Privileges

There was a problem with the privilege set for the ivfServices layout, which made it impossible for the customers to use it. That is now corrected.

PATIENT

Contact

PAT_CONT_ADDITIONALINFO

PAT_CONT_ADDITIONALINFO.KS

New module and variant.

Popover containing additional information about special referral and paper charts.

For more info <http://linnefiler1.foretagssidorna.se/additionalinfo/>



Woman

PAT_WOM_FamilyAnamnesis

PAT_WOM_FamilyAnamnesis.DiseaseAndAge

New module and variant

Possibility to specify family anamnesis for donors.

For more info <http://linnefiler1.foretagssidorna.se/familyanamnesis/>

Partner

Overview sperm analysis

The overview of earlier sperm analyses on the partner page is now customized based on the clinic's choice of module variants on the sperm analysis page. This is also the case for the list of previous sperm analyses on OPU and IUI screen.

TREATMENT

Header information

In layouts related to the treatment table, we have added information in the header area. Around the symbol for treatment, information codes is displayed about the:

- type of stimulation (high level stimulation, natural cycle etc.)
- type of external donation (if it exist)
- planned use or end result of cells (oocyte freeze, thaw of frozen embryos etc.)
- result of the treatment (Cancelled, positive pregnancy, abortion, number of children born etc.).

For more info <http://linnefiler1.foretagssidorna.se/treatment/>

Treatment history

A layout for treatment history is implemented.

Stimulation

Ultrasound performed by and notes at follicle count

The different kinds of stimulation protocols in LinnéFiler have the possibility to switch between a tab for signature of the person who did the measurements of the follicles and a tab for notes. If the user has written a note, the tab will turn red.

Stimulation - thaw

TREATM_ThawStim_DayOfET

TREATM_ThawStim_DayOfET.A

The module variant has changed name to TREATM_ThawStim_CorpLut_at_ETDay.A



Embryo

Fertilisation process & material

ICSI media

Extended popover that gives the possibility to record information about ICSI media.

Embryo freeze

D3

There was a problem that if you froze embryos on day three, the system jumped to the tab of day one. That is now corrected.

Embryo cultivation

Fields that should be hidden were shown in search mode. That is now corrected.

ET

TREATM_ET_ETTypeInfo

TREATM_ET_ETTypeInfo.C

There was a problem with that if you pressed the sum-button the wrong module variant was displayed. That is now corrected.

Problem with a question mark that appeared incorrectly in some fields is corrected.

For more info <http://linnefiler1.foretagssidorna.se/ettypeinfo/>

Sperm analysis

SPA_AfterPrep

SPA_AfterPrep.A

There was a problem with a question mark that appeared incorrectly in some fields. It is now corrected.

For more info <http://linnefiler1.foretagssidorna.se/afterpreparation/>